# Revroll

**Heated Vibrating Foam Roller** 

Model: MG40

**User Manual** 



caring mill by aura

# Warm up. Recover faster. Feel unstoppable.

The Revroll Heated Vibrating Foam Roller blends heat, vibration, and deep compression to melt away muscle knots, boost circulation, and improve flexibility.

From pre-workout prep to post-training recovery — or even injury rehab — it's your personal deep-tissue masseuse, anytime, anywhere.



In the box	5
Unpacking Your Revroll	5 5
	6
Warranty	6
Customer Support	0
Revroll Components	/
Operation Instructions	9
Using Your Revroll	11
Using the Bluetooth Remote Control	12
Vibration Speed Options	13
Using the Heating Pad	13
Tips for Using Your Revroll	15
Maintenance, Safety & Usage Notes	16
Disclaimer	16
Safety Instructions	17
Important Usage Notes	19
Maintenance, Cleaning & Storage	20
Heating Pad Safety Instructions	20

#### In the box

- Vibrating Foam Roller
- Bluetooth Control
- Heated Sleeve
- Battery Charger

#### Unpacking Your Revroll

- Remove your roller from the shipping box and remove all packaging from the unit and accessories.
- Read carefully through this user manual and all safety instructions before using the Revroll and heating pad.
- See Operating Instructions section, in this manual, for information about how to use your Revroll Heated Vibrating Foam Roller.

#### Warranty

The Caring Mill by Aura creates quality, durable health and wellness products and we stand by that quality with a one year warranty on all of our products. If your Revroll ever breaks (and we doubt it will), we'll send you replacement parts and show you how to repair it. If your roller cannot be repaired, we'll replace it, free of charge.

Register your product at www.camar.com/warranty to activate your warranty within 14 days of purchase.

### **Customer Support**

If you have any questions about setting up your Revroll Heated Vibrating Foam Roller, we're here to help. Call or email us and an experienced customer support representative will be happy to assist you: 888.433.6916 or support@camar.com

#### **Revroll Components**

**Heated Sleeve** 

Bluetooth Remote with Wrist Strap & Speed Control



#### **Power and Speed Button**

Press to power the Suger Pro on/off. Then press to toggle through the 5 speed levels

# **Vibration Speed and Battery Level Indicators**

Indicators show current battery level and current vibration speed levels (1-5)

**Charging Port** 



# Special Foam For durability and vibration transfer



# **Diamond Pattern**

The specialized pattern intensifies deep-muscle stimulation



**Charging Adapter** 

#### **Operation Instructions**

Please read the operating instructions before using the Revroll Heated Vibrating Foam Roller. Keep this manual for reference.

#### **Charging the battery**

To charge, connect the AC charging adapter to the indicated port on the foam roller

and to a wall outlet.

Fully charge your roller, for up to 2 hours, before using it for the first time.

After the first charge, your foam roller will take about 3 hours to charge.

Revroll may be recharged at any time, regardless of battery level.



The battery must be charged at least once every 3 months. Do not completely drain the battery before charging.

Average run-time is 2 hours per charge at the highest speed level and 4 hours per charge at the lowest speed level. The power indicator will flash during charging. When the battery is fully charged, the battery indicator will remain lit.

# **USB Output Port**

Plug the heating pad's USB cable into the USB output port to turn it on. The USB Output port can also be used to charge your cellphone or other USB-charging devices

#### **Ventilation Ports**

Do not block ventilation ports.

#### **Using Your Revroll**

- Press the Power/Speed button on the end of the Revroll to power on the device. The Power indicator LED will light up.
- Press the Power/Speed button to toggle through the 5 vibration speed levels: press once for speed level 1, twice for speed level 2, three times for speed level 3, four times for speed level 4, and five times for speed level 5.
- Press the Power/Speed button a sixth time to turn off vibration.
- At the desired vibration speed, roll a muscle group across the The Caring Mill by Aura Foam Roller.
   Apply pressure as needed.
- The unit will auto shutdown after 10 minutes of use. To continue using, just press the Power/Speed button to power on the unit again.

• Always power off the unit after each use.

#### **Using the Bluetooth Remote Control**

- The remote control has an adjustable wrist strap so you can wear it around your wrist for hands-free use.
- Press the Speed button on the remote control five times rapidly to connect it to the Revroll Heated Vibrating Foam Roller via Bluetooth.
- Press the Power/Speed button on the roller to turn it on.
- Press the Speed button on the remote control to toggle through the five vibration speed levels: press once for speed level 1, twice for speed level 2, three times for speed level 3, four times for speed level 4, and five times for speed level 5.
- Press the Power/Speed button on the Revroll to turn off the device.

### **Vibration Speed Options**

Five different massage speeds allow you to choose the right intensity for a particular muscle group or trigger point.

Each speed offers a different vibration frequency (Hz) and percussions per minute (RPM).

Speed	Frequency	RPM
1 2	25 Hz 36 Hz	1,500 RPM 2,160 RPM
3	45 Hz	2,700 RPM
4	53 Hz	3,180 RPM
5	81Hz	4,860 RPM

#### **Using the Heating Pad**

 Connect the heating pad's USB cable to the USB output port on the Revroll.

- Press the Power/Speed button on the Revroll to turn it on and the heating pad will begin to heat.
- Allow the heating pad to pre-heat for about 20-30 seconds.
- Lay on the heating pad to warm up your muscles.
- Begin your massage with the Revroll Heated Vibrating Foam Roller.
- Press the Power/Speed button on the roller to turn it and heating pad off immediately after use.

#### NOTE:

- The heating pad temperature is 122°F.
- The heating pad can only be used with the Revroll. Do not connect the heating pad's USB cable to any other USB output port.

#### **Tips for Using Your Revroll**

- Start with speed level 1 to get your muscles familiar with the roller.
- Slowly roll your body back and forth along
  The Caring Mill by Aura Foam Roller. Take a few
  seconds to go forward and then a few seconds to
  go back. Roll each muscle for several minutes.
- Roll along the long axis of your muscle. For example, if working your triceps, roll from your armpit to your elbow.
- Use a consistent, downward pressure.
- Stay hydrated.
- Use the Caring Mill by Aura Revroll for 15-minute sessions, up to 3 times per day.
- Prolonged application can traumatize tissues.

#### Maintenance, Safety & Usage Notes

#### NOTE:

Please keep this manual in a safe place for reference.

(!) Warning: To reduce risk of burns, fire, electric shock, injury to persons, or damage to equipment, read the following warnings carefully.

#### Disclaimer

The health benefits suggested or implied in this manual, other product literature, and website are not certified or endorsed by any regulatory authority or medical institute.

The use of the vibrating foam roller and heating pad are entirely at the user's discretion. Please read all instructions and safety information carefully before using these products. We assume no responsibility for personal injury or property

damage sustained by or through using these products.

### Safety Instructions

- Warning: For the purposes of recharging the battery, only use the charging adapter provided with this device.
- Consult your doctor before using this product, especially those with diabetes complications, epilepsy, migraines, herniated disks, spondylolisthesis, spondylolysis, or spondylosis, recent joint replacements, and implanted devices such as pacemakers, IUDs, or metal pins or plates.
- Do not use the unit to supplement or replace any medical treatment. Do not use if injured.
- Examine the unit and charging adapter for wear or damage before use. Never use if damaged.

- For adult use only. Keep away from children and pets.
- This device is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision and instruction concerning use of the device by an adult person responsible for their safety.
- Use the device only as instructed in this manual.
   Keep away from liquids and heat sources. Keep ventilation ports free from dust and debris.
- Never use the device while it is charging.
- Unplug the unit after charging and prior to using, cleaning, storing, or transporting. Never leave the device unattended while it is running or charging.
- Do not remove screws or attempt to repair or disassemble the device. The battery is not

removable or replaceable.

 Do not operate the unit continuously for more than one hour. Allow the device to rest for 30 minutes before reusing.

### **Important Usage Notes**

- Use the vibration function only on soft tissues. Do not use vibration on your head, neck, upper back. Do not use on eyes, heart, genitals, or other sensitive areas.
- Check the areas you are treating frequently and stop immediately at the first sign of pain or discomfort.
   Bruising can occur regardless of the vibration speed level selected.
- Use only on dry, clean clothing or skin.
   We recommend using the device over clothing (e.g. form-fitting athletic wear, leggings, etc.). Keep hair, fingers, and other body parts out from under the device.

#### Maintenance, Cleaning & Storage

- Use a slightly damp cloth or soft towel to wipe any dirt or moisture from Caring Mill by Aura Revroll and heating pad after use. Do not use solvents or other chemicals.
- Always power off and unplug the device before cleaning, storing, or transporting the device.
   Store the unit fully charged, in a cool, dry place.
   If storing long-term, charge the battery every
   6 months; holding a charge prevents damage to the battery.

## **Heating Pad Safety Instructions**

 The heating pad may ONLY be used with The Caring Mill by Aura Heated Vibrating Foam Roller. Do not connect the heating pad to any other USB output ports. Keep the USB input dry.

- Turn off the heating pad immediately after use.
   If the heating pad overheats or has a burning smell,
   please stop using it immediately and contact
   customer support.
- Consult your doctor before using this product, especially those with diabetes complications, epilepsy, migraines, herniated disks, spondylolisthesis, spondylolysis,





Revroll
MODEL: MG40
Made in China for Aura Medical LLC